F M A

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510007M01

Introduction

This option is an output of relay contact signal corresponding to judgment results of comparator function of the FMA scale, which can classify into maximum 3 ranks, obtained through D type 25 pin connector. Also, there are signals for “On measurement” and error.

1. Specifications


1-1 pin assignment at connector and the functions (connector: D-SUB25P)

1-1-1 Pins assignment and functions

Pin #	Signal		Input/ Output	Functions • remarks
1	LO	NO	Output	ON when judged as LO
2	LO	NC	Output	OFF when judged as LO
3	LO	COM	-	LO contact common
4	OK	NO	Output	ON when judged as OK
5	OK	NC	Output	OFF when judged as OK
6	OK	COM	-	OK contact common
7	HI	NO	Output	ON when judged as HI
8	HI	NC	Output	OFF when judged as HI
9	HI	COM	-	HI contact common
10	-	-	-	-
11	-	-	-	-
12	-	-	-	-
13	-	-	-	-
14	-	-	-	-
15	-	-	-	-
16	ERR	NO	Output	ON when error is displayed
17	ERR	NC	Output	OFF when error is displayed
18	ERR	COM	-	ERR contact common
19	“On measurement” NO		Output	ON when displayed value is more than 5 times as readability
20	“On measurement” NC		Output	OFF when displayed value is more than 5 times as readability
21	On measurement COM		-	“On measurement” contact common
22	Buzzer	NO	Output	ON when buzzer sounds
23	Buzzer	COM	-	Buzzer contact common
24	External signal (+)		Input	Contact output control(+)
25	External signal (-)		Input	Contact output control(-)

* NO : Normal Open (A contact)

NC : Normal Close (B contact)


D-SUB25P male connector: view from rear

Note


- (1) Be sure to disconnect AC adapter from the power outlet before plugging or unplugging the connector.
- (2) Be sure to use shielded cable for counter measurement against noise.
- (3) To maintain IP65, use a specified water-/dust-proof-type connector. Please contact your local dealer for the specified connector.

1-1-2 Relay contact circuit


WARNING	- When a specified water-/dust-proof-type connector is connected and the scale is to be used in a location where it will get wet during relay contact input, input voltage shall be limited to 16 V~ and 30 VDC .
----------------	---

1-1-3 Remarks on inductive load in the contact


Supply spark-off circuit when inductive load (relay, solenoid etc.) is connected

1-1-4 Equivalent circuit to Input

1-1-5 Example of relay input

Be sure to apply specified current, by avoiding mistake of wrong connection, and using limit resistor and so on. Wrong connection or excessive current supply may damage the internal circuit. Take at least 100 ms for ON period.

2. Setting of the FMA scale

Active the comparator function for outputting the result of comparison through relay contact. Comparator function should be activated even when "ERR" or "On measurement" signal only are to be used.


■ Comparator function setting

Refer to "3-3 Comparator function" in operation manual of FMA.

■ Setting upper and lower limits

Refer to "3-3-2 Setting of the thresholds" in operation manual of FMA.

(example)

3. Operations

1. Plug the AC adapter and press [On/Off] key
2. Enter the set point values (upper & lower limits) to enable the comparator function
(Refer to the previous section and the operation manual of FMA in detail)
3. The Relay contact works as the comparator result is displayed

4. Status of the scale and the relay contact

Status of the scale	Display	Relay contact
Does not meet the discriminant condition	No result is displayed	All OFF
Does not meet the discriminant range	No result is displayed	All OFF
Set points or values are invalid	<LOW>, <OK> and <HIGH> flicker or <ERR748>	All OFF
Over error / Under error	<>>UNDER ERROR<<< / <>>OVER ERROR<<<	- "LO", "OK", "HI" and "On measurement": OFF - "ERR": ON
During setting (not on measurement)	Setting display	All OFF
Waiting for stability with tare or zero	<>>> Tare <<<< or <>>> Zero <<<< flickers	The latest result is held
Stand by	<*>	All OFF
At start up until measurement starts	Start-up display	All OFF

5. Output control by external signal

1. Generally contacts are controlled by judgment result, but can be controlled by external signal by control setting.
2. If power source is connected to external signal input, while output is set to be controlled by external signal ("138 RELAY CTL 2"), judgment result is output for about 200 ms after signal starts.
3. If you apply external signal in pulse wave, take ON period for at least 100 ms.